



ENANTHOL 250

Testosterone Propionate (30mg)

Testosterone Phenyl Propionate (60mg)

Testosterone Isocaproate (60mg)

Testosterone Decanoate (100mg)

FOR INTRAMUSCULAR USE ONLY.

DESCRIPTION:

Sustanol contains 4 different esterifications of natural occurred androgen hormone Testosterone.

As with all the esterifications of testosterone (unless testosterone base) it has a modified release pattern to avoid daily injections.

The testosterone esterifications are used to treat medical conditions such Hypogonadism.

Sustanol 250 is a long lasting complex that allows the administration each 20 days when used for medical conditions treatments.

For fitness purposes the typical procedure is the weekly use or even each 5 days.

This 4 esterifications complex is perfect for medical conditions and fitness uses due to its stable levels of testosterone since the beginning of the application.

ADMINISTRATION:

Sustanol 250 mg is given by deep intramuscular injection only.

DOSAGE:

For hypogonadism treatments the weekly doses are close to 125mg.

For fitness purposes the typical doses in males are 250-500mg weekly.

PRECAUTIONS:

The use of any esterification of testosterone may induce precocious puberty in children and under 20 years patients, if signs of precocious prepuberty occurs the treatment must be discontinued.

Since androgens may cause fluid retention Sustanol 250 must be used with caution in patients with epilepsy, migraine, asthma, and cardiac or renal disease.

ADVERSE REACTIONS:

Headache, irritability, edema, restlessness, depression, tiredness, precocious prepuberty, gynaecomastia, and pain at the site of the injection.

CONTRAINDICATIONS:

Precocious puberty, prostatic carcinoma or another androgen dependent neoplasia, allergy reaction to testosterone or any of the compounds of this formulation.

PRESENTATION:

10ml Multidose vial.

EACH ML CONTAINS:

Testosterone Propionate (30mg).

Testosterone phenyl propionate (60mg).

Testosterone isocaproate (60mg).

Testosterone decanoate (100mg).

Benzyl benzoate , Benzyl alcohol , grape seed oil.

Manufactured for:

NOVATEK LABS

217 Valtaiku iela

Riga - Latvia